Detector de alcool PNI AT6000

- Manual de utilizare -

Atentionari speciale:

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- User manual -

Special warnings:

PNI AT6000 alcohol tester should only be used to identify the presence of alcohol in the breathing air. Do not use this device to determine whether it is safe to drive a vehicle, work with a machine, do dangerous activities or determine the degree of intoxication in a person.

The organism of each person responds differently to alcohol consumption. Therefore, test results with this device should be taken as a reference only.

The manufacturer, importer or retailer does not assume any responsibility for the use of this device. This device should not be used to determine whether a person is capable of driving a motor vehicle in a legal and safe condition. Alcohol consumption, even in small quantities, can affect reflexes and the ability to make decisions.

Before using this device, read the manual carefully and follow the instructions.

Introduction

PNI AT6000 alcohol tester has an advanced sensor with high sensitivity, modern and practical design, being thus a suitable object for personal use. For your safety, when the measured alcohol concentration exceeds the maximum preset level of the device, it will emit a sound and display a warning on the screen.

Main features:

- Alcohol-efficient sensor
- Quick response
- Smart control via MCU (Micro Control Unit)
- Blue illuminated LCD screen
- Display on-screen test process
- Portable design
- Sound warning if the preset limits are exceeded
- Automatic sensor error control
- Low battery indicator

Technical specifications:


- Detection range: 0.00 – 0.20BAC%; 0.00~2.00g/L; 0.00 – 2.00BAC‰~0.00 -1.00 mg/L; 0.00 – 200mg/100mL
- Alarm threshold: 0.05BAC%; 0.50g/L; 0.50 BAC‰; 0.25mg/L; 50mg/100mL
- Accuracy: $\pm 0.02\%$ BAC at 0.05%BAC
- Power voltage: DC4.5V(3×AAA)
- Power consumption: $\leq 120\text{mA}$
- Operating temperature: $-10\text{ }^{\circ}\text{C} \sim 50\text{ }^{\circ}\text{C}$
- Humidity: $\leq 95\%$
- 3-digit illuminated LCD display
- Dimensions: 103 × 65 × 27mm
- Weight: 64 g (without batteries)
- Sensor lifespan: up to 2 ani

Product description and operating mode:

1. Mouthpiece
2. Speaker
3. Button
4. Output for external power
5. LCD screen
6. Ventilation hole

Screen description:

Symbols description:

| | |
|-------|-------------------------------------|
| 8.8.8 | Alcohol content |
| "C" | You may blow |
| "—" | Wait, calculating the concentration |
| FFF | Sensor error |

Basic user instruction:

- Insert into the device 3 AAA 1.5V batteries, respecting polarity.
- Press and hold for one second the on/off button. Immediately, a short acknowledgment sound will be heard. The "Wait" message will be displayed on the screen. Do not blow while the device is heating, otherwise you will get erroneous results.

After the heating period has ended, the device will sound a short acknowledgment sound. The "Blow" message and the "C" symbol will be displayed on the screen. Now you can start testing.

NOTE: The time elapsed before displaying the result depends on the concentration of alcohol. If this is high, the time will be longer.

- Before starting the test, breath air deeply in your chest, then blow directly into the mouthpiece when “C” is displayed on the screen. Stop blowing after you heard a short confirmation sound.
- After you have finished blowing, the symbol “” and the message “TEST” will appear on the screen, which means the device calculates the concentration. Do not blow up the device during measurement. After the concentration analysis is completed, the result will be displayed on the screen. The result will be displayed on the screen for a few seconds, then the device will automatically shut down.

Special warnings:

- Do not cause strong shocks to the device
- If the alcohol concentration in the test air is high, the device may not work normally
- If a test is performed while the battery level is low, the result may not be very accurate
- To ensure the correctness of the test result, wait 15 minutes after you have consumed alcohol
- If the device has not been used for a long time, the result of the first test may be imprecise
- Do not blow directly into the appliance, especially if you are a smoker
- Clean the device periodically. Do not use scouring agents or solvents.

Troubleshooting:

| Defect | Possible cause | Solution |
|--|---------------------------------|---|
| The LCD does not display anything | Battery installed incorrectly | Insert the batteries again, respecting the polarity |
| | Low battery | Replace the batteries |
| | Fired circuits | Contact the seller |
| Lack of gas detection response | The device is not yet warmed up | Wait for the device to warm up completely |
| | Fired circuits | Contact the seller |
| The “FFF” symbol appears on the screen then the device closes | Eroare senzor | Contact the seller |
| The low battery icon appears on the screen, then the device closes | The batteries are very low | Replace the batteries |

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